

Obstetric Labour Table, Economy (Model No. SI/HF/162)



- Overall approximate extended dimension: 1880 mm x 900 mm x 800 mm, 74" x 35" x 31.5" (L x W x H).
- Mattress platform dimension: 1800 mm x 620 mm, 71"x24"(LxW).
- Three section top of S.S. tubular construction.
- Backrest adjustable on ratchet.
- Middle section with "U" cut.
- Leg section can be telescoped under the main frame completely for lithotomy positions.
- Foot operated gas spring assisted trendelenburg position.
- Complete with stainless steel head & side rails hand grips, a pair of rubber padded stainless steel lithotomy rods with rexine ankle straps and stainless steel bowl under 'U' cut.
- Supplied in SKD condition.
- Supplied with suitable rexine covered mattress 25 mm thick & stainless steel I.V. pole.

Optional Accessories

- Lithotomy crutches in lieu of lithotomy rod
- Rexine covered mattress of 25 mm (1") thickness
- Stainless steel I.V. Rod - 2 hooks & 4 hooks