

## Examination-Cum-Gynaec Table (Model No. SI/HF/171)



- Overall approximate dimension: 1830 mm x 610 mm x 825 mm, 72 x 24 x 32.5 (Lx W x H)
- Mild steel tubular construction, machine pressed double bend mild steel sheet
- Two sections with fixed upholstered tops
- Headrest adjusted on ratchet
- Perineal recess and S.S. bowl at leg end
- Pair of synthetic rubber padded stainless steel lithotomy rods with rexine ankle straps
- Legs fixed with rubber feet
- Pre-treated and epoxy powder coated
- Safe working load of 140 kgs and patient load bearing capacity of 135 kgs
- Supplied in knock down condition