

Stretcher Trolley, M.Steel (Model No. SI/HF/149)



- Overall approximate size 2030 mm x 560 mm x 820 mm, 80" x 22" x 32" (L x W x H).
- M.S. tubular framework mounted on four swivel castors,
- 200 mm (8") diameter two with brake and two without brake.
- Removable stretcher top made of aluminium sheet supported on tubular frame.
- Safety belts provided.
- Safe working load of 135 kgs and patient load bearing capacity of 130 kgs.
- Pretreated and powder coated.
- Supplied in knock down condition.

Optional Accessories

- Rexine covered 25 mm thick single mattress 1830 mm x 560 mm, 72" x 22" (L x W)
- I.V. Rod
- Oxygen cylinder cage
- M.S. storage tray