



Blood Pressure Monitor-Arm

Our clients can avail from us a wide array of Blood Pressure Monitor-Arm which is precisely designed using quality tested components and advanced technology. To assure reliable performance and long service life, the offered monitor is stringently tested on predefined quality parameters. Owing to flawless performance, the provided monitor arm is extensively demanded in medical industries. Our esteemed clients can avail this monitor arm from us at market-leading prices.

Features:

- Fuzzy Logic technology
- Arrhythmia detection function (IHB)
- Automatically inflate and deflate
- Automatic power-off function
- Blood pressure classification indicator
- 60 group memory function
- Power: 4 AAA size alkaline batteries
- Measurement Accuracy : Pressure : + 3mmHg Pulse : + 5% of reading