



ANALYTICAL BALANCE

Model : TLP-006

Specification:

General laboratory application modes including totalization, parts counting, display hold...etc.

- With different user selectable weighing units
- Capacity: 200gm
- Readability: 0.01gm
- Repeatability (std.dev.): 0.01g
- Linearity: + 0.02g
- Pan size: 11cm diameter
- Weighing Mode: G, mg, ct
- High contrast display; Auto shut-off
- Power requirement: Battery operated (6 AA batteries) & AC adapter included (electrical)
- Dimension: W x H x D 14cm x 3.5cm x 19cm approximately.