

# GLUCOMETER TEST STRIPS

ACCUCHEK ACTIVE

## ACCU-CHEK

601 B,SILVER UTOPIA, CHAKALA ROAD, ANDHERI (EAST)  
MUMBAI, MAHARASHTRA

Phone: Email: india.accu-chek@roche.com  
<http://www.hospytek.com/accu-chek>

### OVERVIEW:

### DESCRIPTION:

We are widely recognized for providing a wide range of high quality AccuChek Blood Glucose Test Strip. The offered strip is manufactured using the quality assured basic material and latest technique under the guidance of deft professionals. Provided strip is widely used in blood glucose meters to test the blood glucose levels. Moreover, this strip is strictly checked before against numerous parameters so as to avoid any kind of flaws



### TECHNICAL SPECIFICATIONS:

Model Name Active Test Strips - 50 Model Number Active Test Strips - 50 Minimum Blood Sample Needed 1.2 microlitre Measuring Range 10 mg/dL - 600 mg/dL Number of Strips 50 Strips Number of Lancets 0 Lancets Sales Package 1 Vial of 50 Test Strips, 1 Code Chip Other Features Re-dosing Capability, Visual Color Check : Plausibility Check Option with 5 Color Block on the Test Strip Vial ..

READ MORE