

INCUBATOR B.O.D.

18 X 18 X 18

ADAK ELECTRONICS

Village-Unsani Goalbati, Post Office-Unsani, P.S -Jagacha, Nearest Location Mouri Gram Railway Station, Howrah - 711302, West Bengal, India Howrah, West Bengal

Phone: 9804365506 Email: adakelectronics01@gmail.com
<http://www.hospytek.com/adak-electronics>

OVERVIEW:

INCUBATOR

DESCRIPTION:

Double walled construction fully made of St. Steel. The gap between the two walls filled with glass wool insulation. Supplied with S.S. wire mesh shelves. Temperature controlled by Microprocessor based Digital Temp Controller with PT 100 as sensor & safety thermostat. Temp Range : 5°C above ambient to 60°C accuracy $\pm 0.5^\circ\text{C}$



TECHNICAL SPECIFICATIONS:

Forced Air circulating system with fan for better uniformity of Temperature
Electronic Digital
Timer Range 0 to 999 minutes
RS 485 for PC communication with Software Meta
Logger
Single Channel Data logger with printing arrangement ..

READ MORE