

# TWIN OT LIGHT LED

BRS-8

## B.R. SURGICAL

WORK SHOP: 244, STREET NO. 6, NEHRU NAGAR, NEAR ANAND  
PARBHAT INDUSTRIAL AREA ANAND PARBHAT, GALI NO - 10, F  
BLOCK, NEW DELHI - 110008 NEW DELHI, DELHI

Phone: Email: [brsurgical41@gmail.com](mailto:brsurgical41@gmail.com)  
<http://www.hospytek.com/b-r-surgical>

### OVERVIEW:

---

- 4+4 REFLECTOR
- Dome Size:- 500 mm+600 mm
- Main Supply:- 220V /50 Hz

### DESCRIPTION:

---

Controller Fitted with Dome



### TECHNICAL SPECIFICATIONS:

---

|                     |                   |                   |               |
|---------------------|-------------------|-------------------|---------------|
| Illumination:-      | 1,60,000+1,60,000 | Focus:-           | 100 to 200 mm |
| Power Consumption:- | 80+80 W           | No. of LEDs:-     | Single +9X4   |
|                     | 40000-50000 Hrs.  | Validity of Led:- |               |

[READ MORE](#)