

OVERVIEW:

DESCRIPTION:

Perfect Functions It can identify the syringe model automatically (including 10ml, 20ml, 30ml and 50/60ml) and set the speed accurately. It can trigger alarm when occlusion of different levels happens; infusion is near empty and infusion completes. Other functions include battery, mute, rate display, AC/DC power display and so on Convenient Operation It can be easy to single-hand operation by doctors and nurses for unique ergonomics design. The panel of CM-510 is in Chinese and can clearly display about syringe information. Cleanness and Maintenance is convenient. Accurate Infusion The flow rate is accurate and reliable.



TECHNICAL SPECIFICATIONS:

BATTERY back up TOTAL VOLUME DISPLAY;0.1-9999.9ML 10,20,30,50/60ML

SYRINGE Power source : AC: AC220V±10%,50Hz :

Inside battery: chargeable, DC 12V, 1800mAh nickel-hydrogen battery :

Battery time: about 5 hoursSyringe : 10, 20, 30, 50ml syringeFlow rate
range :- 0.1-200ml/h(10ml) 0.1-300ml/h(20,30,50ml)
(0.1ml/h increment)Rapid flow rate :- &n..

READ MORE