

PATIENT MONITOR AMBULATORY / EMERGENCY / EMS MONITOR

C30

COMEN

SHENZHEN,

Phone: Email: info@szcomen.com

http://www.hospytek.com/comen

OVERVIEW:

- COMEN C30 patient monitor with integrated battery
- ECG cable (5 threads) with electrodes
- SpO2 sensor for adults
- Temperature probe
- Power supply
- Standard blood pressure cuff (27-35cm) and hose
- Instructions for use

DESCRIPTION:

The COMEN C30 is a compact, easy-to-operate patient monitor designed to measure ECG, blood pressure, oxygen saturation, heart rate, respiratory rate and body temperature. The shock-proof and waterproof monitor is perfectly suitable for use by rescue service staff and first aid helpers. The COMEN C30 is operated via the touch screen and shortcut keys.

TECHNICAL SPECIFICATIONS:

Compact 4.3" colour touch screen monitor for rescue services
3 modes available: newborn babies, children, adults
ST segments and arrhythmia analysis
48-hour trend memory
With handle
Lightweight and compact
With shortcut keys
Water tight (IPX1 – Protection against dripping water)
Shock-resistant
Battery or mains operation options (approx. 8-hour operating time with battery)
Weight: 1.2kg
Integrated battery
Available with a range of accessories ..

READ MORE