



OVERVIEW:

DESCRIPTION:

The new Compact 525 is now manufactured to deliver a reliable, high quality, more eco-friendly oxygen concentrator.

TECHNICAL SPECIFICATIONS:

Dimensions (H x W x D): 62.2 cm x 34.02 cm x 30.4 cm Weight: 16.3kg Sound level : 40 dbA (50Hz) Flowrate: 0.5 – 5 l / min Oxygen concentration : 93% ± 3% (0.5 – 5 l / min) Electrical requirements : 230V / 50 Hz Power consumption : approx 290 watt at 2 l / min and approx. 312 watt at 5 l / min ..

[READ MORE](#)