

GRAPHIC DISPLAY BIOCHEMISTRY ANALYZERS

CHEM 5X

ERBA (TRANSASIA BIO-MEDICALS LIMITED)

8 Chandivali Studio Road, Andheri (E). Mumbai, MAHARASHTRA

Phone: (022) 40309000 Email: responses@transasia.co.in

<http://www.hospytek.com/erba-transasia-bio-medicals-limited>

OVERVIEW:

DESCRIPTION:

Extensive Bio-chemistry test menu Elaborate Q.C. Program & memory Inbuilt printer Graphic Display 1000 test Result storage



TECHNICAL SPECIFICATIONS:

150 open user-defined test programs
9 Assay modes: 1 point Linear, 2 point Linear, 1 point NonLinear, 2 point Non-linear, Rate A Linear, Rate A Non-Linear, 1 point Sample Blank Linear, 1-point Sample Blank Non-Linear and Absorbance Mono Bichromatic measurement
Photometric range from 0 to 3.0 OD
Photometric resolution up to 0.0001
Static Photometer (340 nm-630 nm)
Wavelength 340 nm, 405 nm, 505 nm, 546 nm, 578 nm, 630 nm
Multi-Point calibration with graphical display
Levy-Jennings graph for quality control
1000 Test Patient Results storage

[READ MORE](#)