

SEMI AUTOMATED BIOCHEMISTRY ANALYZERS

CHEM TOUCH

ERBA (TRANSASIA BIO-MEDICALS LIMITED)

8 Chandivali Studio Road, Andheri (E). Mumbai, MAHARASHTRA

Phone: (022) 40309000 Email: responses@transasia.co.in

<http://www.hospytek.com/erba-transasia-bio-medicals-limited>

OVERVIEW:

DESCRIPTION:

Real Time Graph Touch Screen Interface
Low Reagent Consumption
Triple Cuvette System
USB Conn/activity to Printer PC
Built in Self Diagnosis



TECHNICAL SPECIFICATIONS:

Static Photometer with photometric range (340 nm to 670 nm.)
8 Interference Filter 340 nm, 405 nm, 505 nm, 546 nm, 578 nm, 600 nm, 670 nm.
Silicon Photo Diode detectors
Photometric Range from 0 to 3.0 O.D.
Quartz Halogen Lamp of 12V, 20W
High-Resolution Graphics LCD with Backlight
View Area: 120x92 mm ..

[READ MORE](#)