

INT KNEE O DRAPE

I525

G.Surgiwear Limited

Rasoolpur Jahanganj Near Hathoda Crossing, Shahjahanpur - 242001
India Shahjahanpur, Uttar Pradesh

Phone: +91 5842223818 Email: mukherjee_surgiwear@yahoo.co.in,
GSL_Lucknow@surgiwear.in
<http://www.hospytek.com/g-surgiwear-limited>

OVERVIEW:

DESCRIPTION:

SURGIWEAR has developed this drape after feedback from surgeons. It offers many advantages: Structured construction of this drape provides: Level 4 protection as per ANSI / AAMI Standards. No fluid can pass through the drape in operating area. Minimum transfer of fabric fibres into the wound. No chemicals present in the fabric of the drape.



TECHNICAL SPECIFICATIONS:

For all kind of surgery in knee & below area including knee replacement. Material: 65-77 GSM Spun lace fabric water repellent Size: 240x310 cm, Elastic aperture: 5 cm dia ..

[READ MORE](#)