

INT \ ELBOW O DRAPE

I526

G.Surgiwear Limited

Rasoolpur Jahanganj Near Hathoda Crossing, Shahjahanpur - 242001
India Shahjahanpur, Uttar Pradesh

Phone: +91 5842223818 Email: mukherjee_surgiwear@yahoo.co.in,
GSL_Lucknow@surgiwear.in
<http://www.hospytek.com/g-surgiwear-limited>

OVERVIEW:

DESCRIPTION:

SURGIWEAR has developed this drape after feedback from surgeons. It offers many advantages: Structured construction of this drape provides: Level 4 protection as per ANSI / AAMI Standards. No fluid can pass through the drape in operating area. Minimum transfer of fabric fibres into the wound. No chemicals present in the fabric of the drape.



TECHNICAL SPECIFICATIONS:

For all kind of surgery on elbow and forearm
Material: 65-77 GSM Spun lace fabric water repellent
Overall size: 240x220 cm, Elastic aperture: 5cm dia, Poly sheet lining: 120x120 cm

[READ MORE](#)