

# SHOULDER ARTHROSCOPY DRAPE

D509

**G.Surgiwear Limited**

Rasoolpur Jahanganj Near Hathoda Crossing, Shahjahanpur - 242001  
India Shahjahanpur, Uttar Pradesh

**Phone:** +91 5842223818 **Email:** mukherjee\_surgiwear@yahoo.co.in,  
GSL\_Lucknow@surgiwear.in  
<http://www.hospytek.com/g-surgiwear-limited>

## OVERVIEW:

---

## DESCRIPTION:

---

SURGIWEAR has developed this drape after feedback from surgeons. It offers many advantages: Structured construction of this drape provides: Level 4 protection as per ANSI / AAMI Standards. No fluid can pass through the drape in operating area. Minimum transfer of fabric fibers into the wound. No chemicals present in the fabric of the drape.



## TECHNICAL SPECIFICATIONS:

---

Material: 65-77 GSM Spun lace fabric water repellent  
Size: 160x190 cm(1 sheet), 160x210 cm(1 sheet)..

[READ MORE](#)