

BP MONITOR BP CUFF

JMS-BP

JHALANI MEDICAL & SURGICAL

E-66, Dawa Bazaar, Madhav Club Road, Freeganj, Ujjain - 456010
Freeganj, Madhav Nagar

Phone: 0734-2550167 Email:

<http://www.hospytek.com/jhalani-medical-surgical>

OVERVIEW:

DESCRIPTION:

Home BP should be measured using a validated, automatic BP device (preferably with memory storage), using an appropriately sized upper arm cuff. Measurements should be taken after five minutes of seated rest and before medication, food or vigorous exercise. BP should be recorded for seven days (five days minimum) in the morning and evening (two readings each). Overall, home BP is the average systolic and diastolic BP over seven days (excluding the first day); an average of $\geq 135/85$ mmHg is indicative of hypertension.



TECHNICAL SPECIFICATIONS:

Type-Upper Arm Battery Type-AA Number of Batteries Required-4 Minimum Pressure Measurement Range-0 mmHg Maximum Pressure Measurement Range-299 mmHg Minimum Pulse Measurement Range-40 beats/min Maximum Pulse Measurement Range-180 beats/min Pulse Measurement Accuracy: ± 5 Pressure Measurement Accuracy: ± 3 Sales Package-Upper Arm BP Monitor Device ..

READ MORE