

PORTABLE NERVE MUSCLE STIMULATOR

STIMTEC 4

JOHARI DIGITAL HEALTHCARE LTD

G-582-584, EPIP, BORANADA- SALAWAS ROAD, BASNI SILAWATAN, RAJASTHAN Basni Silawatan, Rajasthan

Phone: +91-291-3091104 Email: info@joharidigital.com

<http://www.hospYTEK.com/johari-digital-healthcare-ltd>

OVERVIEW:

- Symptomatic relief of chronic and intractable pain.
- Adjunctive treatment of post-surgical and post-traumatic acute pain conditions.
- Portable handheld electrotherapy device.
- Maintaining or increasing range of motions.

DESCRIPTION:

Stimtec 4 is a portable electrotherapy device featuring 4 therapeutic modes: TENS, Interferential, EMS, Russian. This handy device is used for pain relief and electrical muscle stimulation.



TECHNICAL SPECIFICATIONS:

Main power supply	230 Vac, 50-60Hz, $\pm 10\%$ - 115 Vac, 50-60 Hz*, $\pm 10\%$ Channels
	2WaveformsTENS/EMS
	Mono-phase square pulse. IF/RUS
	Synchronous, Alternate, Delay Russian ModesIF &nb ..

[READ MORE](#)