

ULTRASOUND CAVITATION

HI-LIPO CAVI

JOHARI DIGITAL HEALTHCARE LTD

G-582-584, EPIP, BORANADA- SALAWAS ROAD, BASNI SILAWATAN,
RAJASTHAN Basni Silawatan, Rajasthan

Phone: +91-291-3091104 **Email:** info@joharidigital.com

<http://www.hospytek.com/johari-digital-healthcare-ltd>

OVERVIEW:

DESCRIPTION:

Hi-Lipo Cavi is the fast and non-surgical type of lipo-suction device. It uses low frequency 40KHz Ultrasound treatment for non surgical reduction of localized fat. The body areas treated by Hi Lipo Cavi react in depth with the formation of bubbles that implodes, causing the rupture of adipocytes, which will then be disposed off by the liver cells and renal besides a good physical activity.



TECHNICAL SPECIFICATIONS:

Main power supply - 230 Vac, 50-60Hz, $\pm 10\%$ - 115 Vac, 50-60 Hz*, $\pm 10\%$ Output - 1
(Ultrasound Cavitation)Modes - Continuous and PulseFrequency - 40 KHzPower Output - 50
wattsTimer - 30, 45 & Continuous..

[READ MORE](#)