

OVERVIEW:

DESCRIPTION:

JUS-2 provides the most effective Continuous and Pulsed treatments in pain relief and muscle rehabilitation . It delivers 1 & 3 MHz through a unique ultrasound applicator. The effect of ultrasound via an increase in local blood flow can be used to help reduce local swelling, chronic inflammation and promote bone fracture healing.



TECHNICAL SPECIFICATIONS:

Main power supply: 110-220 V AC Mode: Pulsed, Continuous Frequency: 1 MHz & 3 MHz % Duty Cycle: 10, 20 & 50% Output: One Max. Intensity: 3 watt/cm² Treatment Time: 1 – 60 Minutes Electrical Type: TYPE B Shipping Weights: 4.5 Kg Shipping Dimension: 41 x 38 x 20 cm (LWH) ..

[READ MORE](#)