

COMBINATION THERAPY UNIT

WINSTIM WS2U

JOHARI DIGITAL HEALTHCARE LTD

G-582-584, EPIP, BORANADA- SALAWAS ROAD, BASNI SILAWATAN,
RAJASTHAN Basni Silawatan, Rajasthan

Phone: +91-291-3091104 Email: info@joharidigital.com

<http://www.hospytek.com/johari-digital-healthcare-ltd>

OVERVIEW:

DESCRIPTION:

WINSTIM is a clinical combination therapy of Ultrasound, Electrotherapy with easy user interface and versatility to treat different body areas simultaneously. This aesthetically designed clinical model has a unique ultrasound head with inbuilt power control and user preset electrotherapy programs. The state of the art WINSTIM is light weight, small in size, and battery powered.



TECHNICAL SPECIFICATIONS:

Up to 50 Pre-programmed protocols Clinical Combination Therapy. Electro diagnostics: I/t curves, rheobasis and chronaxie Full range of low and medium frequency waveforms Relaxation of muscle spasm and Joint contracture. Simultaneous treatment of multiple patients. Ultrasound under Water Therapy for Head Warming. Capacitive touch screen LCD. US FDA 510(k) cleared. 45-minute battery back-up. Advanced Smart-Phone style OS..

READ MORE