

# SURGICAL LIGHT LED

XN-233

## MEDITECH INDIA

487/67, Glass Street, Opposite Punjab National Bank, National Market,  
Peeraghari New Delhi - 110087, Delhi, India New Delhi, Delhi

Phone: (11) - 25255415 Email: [marketinguk@meditech-india.com](mailto:marketinguk@meditech-india.com)  
<http://www.hospytek.com/meditech-india>

### OVERVIEW:

### DESCRIPTION:

We offer 2nd Generation LED Surgical light, innovative thought to build up operating room with powerfull illumination light for better surgical. Our LED Surgical light provides cold, shadowless and brighter light area for long surgeries, without any interuption to surgeons. LEDs have a technical advantage over the traditional halogen lamp, which produce heat, limited life cycle of bulb. Thanks to LEDs technique, now surgeon can work in low heat and great life cycle of LEDs working hours, with low power consumption & natural white day light



### TECHNICAL SPECIFICATIONS:

Model No. : XN 233  
Light head dia: 53 cm x 2  
Illumination intensity @1 m: 250,000 Lux  
Focus dia : 15-20 cm  
Color temperature: ..

[READ MORE](#)