

TROLLEY BABY BASSINET

50-1000 D

Narula Udyog India Private Limited

NARULA'S TOWER - I, 25, CENTRAL MARKET, PUNJABI BAGH
(WEST), NEW DELHI – 110026, INDIA New Delhi, Delhi

Phone: +91-11-2522 9526, 4246 3777 **Email:** sales@medikraft.com
<http://www.hospytek.com/narula-udyog-india-private-limited>

OVERVIEW:

- Baby Bassinet on Trolley

DESCRIPTION:

S.S. Tubular Framework mounted on 10 cms Castors (Two with Brakes).
Transparent Bassinet for the infants. Provided with: 3 cms thick Mattress, I.V.
Rod, Mosquito Net Holder and a Utility Tray at the base.



TECHNICAL SPECIFICATIONS:

[READ MORE](#)