

PULSE OXIMETER TABLE TOP

NL 5100

NECTAR LIFESCIENCES LTD

SCO 38-39, SECTOR 9 D. Chandigarh, Punjab

Phone: 172-3078158 Email: marketing@neclife.com

<http://www.hospytek.com/nectar-lifesciences-ltd>

OVERVIEW:

DESCRIPTION:

Routine Check mode & Continuous check mode
Data graph & Trend Table review
Perfect Table-mount solution
High SD Memory



TECHNICAL SPECIFICATIONS:

Display 2.8 (320-240) color TFT LCD and red and yellow highlight LED. Visual and audible alarm for Systolic Pressure(SYS), Diastolic Pressure(DIA), Mean Pressure(MAP), SpO2 and Pulse rate(PR) Unique Flash memory, storage for up to 2000 NIBP data and 78000 SpO2 data..

[READ MORE](#)