

WHEEL CHAIR ELECTRIC

JAZZY J6

NUK HEALTHCARE

Shop No. 12, Sonal Apartments, Bldg. No.2 Sonawala Cross Road
No.1, Opp.SBI, Goregaon (East), Mumba Mumbai, Maharashtra

Phone: 9820207899 Email: nitin@nukindia.com
<http://www.hospytek.com/nuk-healthcare>

OVERVIEW:

DESCRIPTION:

Drive Wheels - 10" Solid 110" Pneumatic2 Caster Wheels Front: 5" Solid Rear: 6" Solid Maximum Speed ^4/^ Up to 4 mph Ground Clearance 2.25" (battery trail rail) Turning Radius ^5/^ 21.25" Overall Length 34.5" without foot rest Base Width 23.25" Seating Sizes: Width / Depth 10" - 20" Seat-to-Floor Heights: 16.5" - 17.5" Seat-to-Floor Heights: 16.875" Standard Battery Size ^6/^ U-1 Battery Weight ^14/^ 24.5 lbs. Standard Charger 5A Off-board Weight Capacity ^7/^ 300 lbs Base Weight 104.5 lbs. Available Electronics ^12/^ 70A Q-Logic NE 70A Q-Logic NE+ 70A Q-Logic EX



TECHNICAL SPECIFICATIONS:

READ MORE