

STRETCHER DOUBLE FOLD

PH-134

PRESENS

India Road,


Phone: Email:

<http://www.hospytek.com/presens>

OVERVIEW:

DESCRIPTION:

Stretcher Double Fold, Light Weight, easy clean, anodized aluminium alloy framework. Washable, water proof & anti-fungus fabric top with safety belts,



TECHNICAL SPECIFICATIONS:

Unfolded dimensions (L x W x H): 213 x 56 x 15 cm
Folded dimensions (L x W x H): 110 x 20 x 11 cm
Self weight: 7.5 kg
Load bearing: 140 kg

[READ MORE](#)