

OXYGEN CONCENTRATOR DOUBLE FLOW

OXYCURE 240 A WITH NEBULISER

S-CURE

Phone: Email:

<http://www.hospytek.com/s-cure>

OVERVIEW:

DESCRIPTION:

Easy maintenance
Optional wifi- remote data transport
Low oxygen purity alarm
Continuous delivery method



TECHNICAL SPECIFICATIONS:

Oxygen purity 93 ±3%
O₂ purity alarm : low oxygen at 82%
Metal connector for fire proof
Flow Rate Oxygen flow: 0.5 -5 L/min ..

[READ MORE](#)