

HOSPITAL BED SEMI FOWLER

MS1008

SURGICAL DEVICES HUB

No. 26, Netaji Subhash Marg, Near Mc Donalds Daryaganj, New Delhi - 110002, Delhi, India Delhi, DELHI

Phone: Email: mahajanexport@gmail.com
<http://www.hospytek.com/surgical-devices-hub>

OVERVIEW:

HOSPITAL BED SEMI FOWLER have special features both for the comfort and well being of patient and convenience of healthcare workers. Hospital Bed is specially designed for hospitalized patients or in a need of some type of healthcare .

DESCRIPTION:

Mechanically operated, uniformly perforated two sectional press bent CRCA seat top bed. Preciserectangular steel tubular frame work. One fixed levelhandles at foot end to obtain positions (like head raise)by smooth crank mechanism. Detachable head and footABS (IMPORTED) Bows. With IV stand epoxy powdercoated.



TECHNICAL SPECIFICATIONS:

Model No. : MS1008Type of Product : BedSub Type : Patient Semi fowler ..

READ MORE