

PREMIER LATEX NATURAL RUBBER GLOVES

IND-962- 6.5

SURGICARE

50A, Senguptha Street, Ramnagar, Coimbatore - 641009, Tamil Nadu,
India Coimbatore, Tamil Nadu

Phone: 0422-2230023 , 0422-2232223 Email:

surgicare.cbe@gmail.com

<http://www.hospytek.com/surgicare>

OVERVIEW:

DESCRIPTION:

Surgicare Premier is made from Natural Rubber Latex. The glove is anatomically shaped with thumb ball effect, reducing hand fatigue during surgical procedures. Polymercoated Surgicare provides excellent barrier protection for both doctor and patient.



TECHNICAL SPECIFICATIONS:

TYPE: Powder Free Polymer Coated Latex Surgical Gloves. Sterilized by ETO.PRIMARY
MATERIAL: Natural Rubber Latex.DESIGN: Hand Specific, Curved finger, Beaded cuff,
Micro rough textured.POWDER FREE RESIDUE: < 2 mg/gloveSTERILIZATION: Ethylene
Oxide / Gamma as per Customer requirement..

[READ MORE](#)